




Prevantics® Swab2Hub™ is a disinfectant for the external parts of needleless connectors, hubs and ports.

- + **Complies with Guidelines:** its 2% Chlorhexidine gluconate (CHG) 70% Isopropyl alcohol (IPA) formulation is recommended in evidence-based infection guidelines published by health organisations such as NICE<sup>1</sup>, epic3<sup>2</sup>, and CDC<sup>3</sup>.
- + **Proven Efficacy:** delivers a broad spectrum of antimicrobial action with excellent residual activity proven to reduce catheter related bloodstream infections by 75%<sup>4</sup>.
- + **Consistent & reliable results:** reduced guess work, providing clinicians with a point of use solution, and standardised best practise for optimal infection prevention.
- + **Designed for convenience:** Ready to use, single-use and easy to open.



Designed for consistency:

By implementing protocols and standardising practises for the disinfection of needleless access sites, your facility can help reduce the risk of contamination and improve outcomes.



Insert the needleless connector end into the cavity of the Prevantics® Swab2Hub device and twist in circular movements for at least 30 seconds.

The innovative designed device provides friction to scrub and disinfect, removing microorganisms from the grooves and threads of needleless connectors.

Product efficacy:

SWAB2HUB	TEST METHOD	ORGANISM	CONTACT TIME*	
Bacteria	EN 13727	Pseudomonas aeruginosa	30 sec (c)	
	EN 13727	Staphylococcus aureus	30 sec (c)	
	EN 13727	Enterococcus hirae	30 sec (c)	
Yeasts	EN 13624	Candida albicans	30 sec (c)	
Enveloped Virus	EN 14476	Vaccinia Virus	15 sec (c)	60 sec (d)
	EN 14476	Corona,HIV, Hep B & C, SARS, Flu	15 sec (c)	60 sec (d)
TB	EN 14348	Mycobacterium terrae	60 sec (c)	

\*Clean Conditions (c) / Dirty Conditions (d)

	REORDER NO	QUANTITY PER CARTON	QUANTITY PER CASE	CLASSIFICATION
Swab2Hub	ISH4000	1x100	8 (800 units)	Medical Device Class IIb



Ready to use disposable disinfection wipes for medical devices containing Chlorhexidine gluconate (CHG) at 2% and Isopropyl alcohol (IPA) at 70%.

## Product efficacy:

2% CHG SWAB	TEST METHOD	ORGANISM	CONTACT TIME*	
<b>Bacteria</b>	EN 13727	Pseudomonas aeruginosa	30 sec (c)	
	EN 13727	Staphylococcus aureus	30 sec (c)	
	EN 13727	Enterococcus hirae	30 sec (c)	
<b>Yeasts</b>	EN 13624	Candida albicans	30 sec (c)	
<b>Enveloped Virus</b>	EN 14476	Vaccinia Virus	15 sec (c)	60 sec (d)
	EN 14476	Corona, HIV, Hep B & C, SARS, Flu	15 sec (c)	60 sec (d)
<b>TB</b>	EN 14348	Mycobacterium terrae	60 sec (c)	

\*Clean Conditions (c) / Dirty Conditions (d)

## Range details:

	REORDER NO	WIPE COUNT	WIPE SIZE	QTY PER CASE	CLASSIFICATION
200ct Sachet	ISW4050	200	150mm x 129mm	8 (1600 units)	Medical Device Class IIa
100ct Sachet	ISW4040	100	150mm x 129mm	12 (1200 units)	Medical Device Class IIa

## Other formulations available:



## References:

<sup>1</sup> National Institute for Clinical Excellence. Clinical Guideline 2. Infection Control. Prevention of healthcare-associated infection in primary and community care. June 2003. <http://www.nice.org.uk/pdf/CG2fullguidelineinfectioncontrol.pdf>

<sup>2</sup> Loveday, HP, et al., 2014. Epic3: National Evidence-Based Guidelines for Preventing Healthcare-Associated Infections in NHS

Hospitals in England. J. Hosp. Infect. 86:S1-S70.

<sup>3</sup> Centers for Disease Control and Prevention. (2002) Guidelines for the prevention of intravascular-catheter-related infections. MMWR 51 (RR-10): 1-29. [www.cdc.gov/mmwr/PDF/rr/rr5110.pdf](http://www.cdc.gov/mmwr/PDF/rr/rr5110.pdf)

<sup>4</sup> Soothill JS, Bravery K, Ho A, Macqueen S, Collins J, Lock P. "A fall in bloodstream infections followed a change to 2% chlorhexidine in 70%

isopropanol for catheter connection antisepsis: A pediatric single centre before/after study on a hemopoietic stem cell transplant ward". Am J Control 2009 Oct;37-8:626-30.

<sup>5</sup> Disinfection of Needleless Connector Hubs: Clinical Evidence Systematic Review. Nancy L. Moureau and Julie Flynn. Nurs Res Pract. 2015; 2015: 796762.

**Find out more about how our Layered Approach to Infection Prevention can help you.**

PDI EMEA, Pywell Road, Willowbrook Est. Industrial Estate, Corby - Northants, NN17 5XJ, UK

+44 (0) 8081 697 945   contact@pdi-emea.com   wearepdi-intl.com   PDI International   PDI\_Intl

