













Super Sani-Cloth® is a two in one solution ideal for cleaning and disinfecting in fast-paced environments that require short contact times and broad microorganism efficacy.

- + **Ready to use:** pre-dosed, premium, strong, low linting disposable cleaning and disinfection wipe, requiring no mixing or measuring, with no smear formulation.
- + **Proven**[efficacy: effective against Coronavirus and Norovirus in as little as 60 seconds.
- + **Improved Efficiencies:** one product to clean and disinfect, providing you with consistent and reliable results.
- + **Superior Compatibility:** a quaternary ammonium, alcohol-based formulation which can be used on non-invasive medical devices; non-porous hard surfaces with supporting compatibility data available.
- + **Cost** effective: using the same solution for cleaning and disinfecting saves you time and money.





Product Efficacy

SUPER SANI-CLOTH®	TEST METHOD	ORGANISM	CONTACT TIME*
Bacteria	EN 16615	Pseudomonas aeruginosa	60 sec (c) 60 sec (d)
	EN 16615	Staphylococcus aureus	60 sec (c) 60 sec (d)
	EN 16615	Enterococcus hirae	60 sec (c) 60 sec (d)
	EN 13727	Pseudomonas aeruginosa	30 sec (c)
	EN 13727	Staphylococcus aureus	30 sec (c)
	EN 13727	Enterococcus hirae	30 sec (c)
Yeasts	EN 16615	Candida albicans	60 sec (c) 60 sec (d)
	EN 13624	Candida albicans	30 sec (c)
	EN 13624	Candida auris	30 sec (c)
Enveloped Virus	EN 14476	Vaccinia Virus	60 sec (c)
	EN 14476	Corona,HIV, Hep B & C, SARS, Flu	60 sec (c)
Non-enveloped Virus	EN 14476	Norovirus	60 sec (c)
ТВ	EN 14348	Mycobacterium terrae	60 sec (c)

*Clean Conditions (c) Dirty Conditions (d)

Point of Use Accessories



Bracket for canisters Bracket: A001PUR

Canister brackets affix securely to the wall, holding the canister in place at its location.



Pole Adapter Product Code: A002POL

Pole brackets make it possible to affix the canister brackets on to portable equipment such as BP stands.



I Am Clean Sticky Notes Product Code: A005SNT

Recyclable sticky notes which can be used to easily and effectively identify equipment or surfaces that have been



Educational wall charts and signage

Product evaluations, audits, educational wall charts & signage, that can be personalised for your specific

Range Details

	REORDER NO	WIPE COUNT	WIPE SIZE	QTY PER CASE	CLASSIFICATION
150ct Short Canister	ESC1630	150	130mm x 130mm	6	Medical Device Class IIa
100ct Large Canister	ESC1620	100	200mm x 195mm	6	Medical Device Class IIa

Other formulations available:







Find out more about how our Layered Approach to Infection Prevention can help you.

PDI EMEA, Pywell Road, Willowbrook Est. Industrial Estate, Corby - Northants, NN17 5XJ, UK











