



# General Guidelines For Use

## PT1 Human Hygiene Biocide.

Please follow local healthcare policy for personal protective equipment, hand hygiene & waste disposal guidelines.



**1.** Safely open the packaging without any external contact with the product.



**2.** Safely remove the swabstick from the packaging without any external contact, assisting aseptic non-touch technique (ANTT) compliance.



**3.** The 2-sided applicator can be used around joints, interdigital areas and skin folds.



**4.** Apply the liquid to the desired area of the skin to be disinfected. Allow to air dry.

**Find out more about how our Layered Approach to Infection Prevention can help you.**

PDI EMEA, Pywell Road, Willowbrook Est.  
Industrial Estate, Corby - Northants, NN17 5XJ, UK

+44 (0) 8081 697 945

contact@pdi-emea.com

wearepdi-intl.com

PDI International

PDI\_Intl

Space for hospital logo